


FIEBRE Standard Operating Procedure**Title** Performing the Direct Agglutination Test (DAT) for Leishmania**SOP Reference****Version****Date of effect**

1.1.1

SOP Development

	Name	Title	Signature	Date
Author	Victoria Gould	LSHTM Laboratory Co-ordinator	Victoria Gould	02/03/21
Reviewer	Chrissy h Roberts	LSHTM FIEBRE laboratory expert		02/03/21
Approver	Heidi Hopkins	FIEBRE scientific program coordinator	Signed by email	17/02/21

Review Tracker

Due date for next review	Reviewer name	Signature	Date reviewed

Revision History

Version No.	Effective date	Reason for change

SOP User Confirmation

I acknowledge that I have read, understood and agree to follow this SOP

#	Name (print)	Signature	Date
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			

1. Title: Performance of the Direct Agglutination Test (DAT) for Leishmania

2. Purpose: To describe the procedure for carrying out the DAT on serum samples from the FIEBRE study. This SOP describes procedures for selection, dilution and assaying of serum samples for agglutination with the DAT assay to determine positivity for Leishmaniasis infection.

3. Responsible staff: FIEBRE laboratory staff at LSHTM

4. Background & Rationale: The first primary objective of the FIEBRE study is to determine the treatable and/or preventable causes of fever in our study population. One of the potential causes in the study areas are Leishmania species. The DAT assay is well established and works on the principle that where antibodies are present agglutination with Leishmania antigen occurs providing a visible readout.

5. Supplies and Materials

- Serum sample aliquot allocated for Leishmania testing
- Freeze dried, Coomassie blue stained, Leishmania antigen
- V bottom microtitre plates
- Saline
- Urea saline
- Beta Mercaptoethanol
- Reagent reservoirs
- Pipettes (including Multichannel)
- Suitable pipette tips

- Biohazard disposal system
- Lab surface disinfectant (e.g. Virkon)
- Sterile filtration units
- Heat block

6. Procedures:

6.1 Sterilisation of solutions

6.1.1 Make saline to 1x solution (if required).

6.1.2 Filter through a unit containing 0.2µM pore.

6.2 Reconstitution of *Leishmania antigen*

6.2.1. Reconstitute the freeze-dried antigen in 5mL of sterile saline and rotate (**DO NOT SHAKE**) to ensure a homogenous solution

6.3 Selection of samples

6.3.1 Carry out initial testing on day 28 samples where there is sufficient sample to allocate for Leishmania testing. If positive, or where no sample is available, test the day 0 sample.

6.3.2 Remove the test samples from storage at -80°C and thaw prior to testing. Test the samples in the Cryptococcal antigen test, if required, then heat inactivate the sample at 56°C. Place samples in the fridge after use to facilitate the titration of any positive samples as soon as possible without any additional freeze thaw cycles.


6.4 Dilution of samples

6.4.1 Mix a suitable quantity of urea saline with Beta Mercaptoethanol (70 µL of 2-ME per 10 mL of saline) and pour into a reagent reservoir for use as sample diluent.

6.4.2 Initially a single sample dilution will be screened. Carry out a Pre-dilution step by adding 318µL of diluent to all wells of a dilution plate using a multichannel pipette and suitable tips.


6.4.3 Add 2 µl of sample to wells in columns 2-12 and mix. Add positive controls to each plate in appropriate wells in column 1, along with a negative control. Record the ID of samples added to each well in the test template (Appendix A). *See plate map.*

6.4.4 Add 45µl of sample diluent to the test plate. Transfer 5 µl from the dilution plate to corresponding wells in the test plate and mix. This plate now contains the screening dilution of 1:1600.



6.4.5 For positive samples at the screening dilution carry out a sample titration by adding 98 μ L of diluent into wells A1 to H1 and 50 μ L into the remaining wells. Add 2 μ L of sample to wells in column 1 along with controls.

6.4.6 Remove 50 μ L from column 1 and serially dilute through columns 2-11. Discard the extra 50 μ L from column 11. Column 12 is left blank. See *plate map*




7 Addition of Antigen

- 7.1.1 Add 50 μ L of antigen to all wells. Start with the blank wells to minimise the risk of cross contamination.
- 7.1.2 Incubate the plate for a minimum of 18 hours at room temperature, undisturbed, for any antibody/antigen reaction to occur

8 Reading of results

- 8.1.1 After the incubation period examine the plate for signs of agglutination, where there will be a uniform blue colour across the well in comparison to a negative well where the unreacted promastigotes fall to produce a blue button in the point of the well.
- 8.1.2 Where titrations are being done the titre is the last dilution that shows a difference in agglutination compared to the blank or negative control.
- 8.1.3 There is a cut off of 1:1600. Report titres of 1:100 to 1:800 as negative. Report titres of 1:1600 or above as positive.
- 8.1.4 Example plate (NB antigen control is column 1 not 12 in this example)



9 Reporting of results

- 9.1.1 Enter the results on to the DAT test template
- 9.1.2 Report the results using the appropriate form (e.g on REDCap)

10 Appendix

Leishmania DAT test template

Date _____ Set Up By _____ Ag Lot No. _____

Urea Saline Lot No. _____ Positive Control _____ Negative Control _____

Read By/Date _____ Check by/Date _____

Plate No	1:100	1:200	1:400	1:800	1:1600	1:3200	1:6400	1:12800	1:2560	1:512000	1:102400	0
	1	2	3	4	5	6	7	8	9	10	11	12
Positive Control												
Negative Control												
Sample 1												
Sample 2												
Sample 3												
Sample 4												
Sample 5												
Sample 6												

Plate No.	1	2	3	4	5	6	7	8	9	10	11	12
A	+ve Ctrl	1	9	17	25	33	41	49	57	65	73	81
B	+ve Ctrl	2	10	18	26	34	42	50	58	66	74	82
C	Low +ve	3	11	19	27	35	43	51	59	67	75	83
D	Low +ve	4	12	20	28	36	44	52	60	68	76	84
E	-ve Ctrl	5	13	21	29	37	45	53	61	69	77	85
F	-ve Ctrl	6	14	22	30	38	46	54	62	70	78	86
G	blank	7	15	23	31	39	47	55	63	71	79	87
H	blank	8	16	24	32	40	48	56	64	72	80	88

Sample	ID	Result	Comment
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			
20			
21			
22			
23			
24			
25			
26			
27			
28			
29			
30			
31			
32			
33			
34			
35			
36			
37			
38			
39			
40			
41			
42			

43			
44			
45			
46			
47			
48			
49			
50			
51			
52			
53			
54			
55			
56			
57			
58			
59			
60			
61			
62			
63			
64			
65			
66			
67			
68			
69			
70			
71			
72			
73			
74			
75			
76			
77			
78			
79			
80			
81			
82			
83			
84			
85			
86			
87			
88			